



Specialised
Orthotic
Services

MoJo & Stroller 2 Net Carrier



USER'S MANUAL

ALL USER'S OF THE EQUIPMENT SHOULD BE AWARE
OF THIS DOCUMENT AND ITS CONTENT

CONTENTS



1 x Net Carrier
1 x Base Board (pre-installed)

Note:

If you remove the base board for cleaning, etc make sure that the narrow end is fitted to the rear of the net carrier (label end) when refitting.

1. MOJO NET CARRIER FITTING INSTRUCTIONS

- Place the bag inside of the chassis with the manufacturers label at the rear by the tilt pedal.
- Hook the rear straps over the rear cross beam as shown (Fig.1a) and secure the two halves of the Velcro together.
- Hook the front straps over the front cross beam as shown (Fig.1b) and secure the two halves of the Velcro together.
- Hook the middle straps over the lower frame, around the back of the tilt mechanism as shown (Fig.1c) and secure the two halves of the Velcro together.
- Correctly fitted net carrier shown (Fig.1d).
- Removal is a reversal of the above.
- **Be careful not to overload the carrier especially with loose items as these may lodge in the wheels causing instability.**



Fig.1a



Fig.1b



2. STROLLER 2 NET CARRIER FITTING INSTRUCTIONS

- Place the carrier inside the chassis with the manufacturers label at the rear by the tilt pedal.
- Hook the rear straps over the rear cross beam as shown (Fig.2a) and secure the press studs together.
- Hook the side straps over the frame as shown (Fig.2b & 2c) and secure the two halves of the Velcro together.
- Correctly fitted net carrier shown (Fig.2.1c).
- The rear section of the net carrier can be expanded via the two zips.
- Removal is a reversal of the fitting procedure.
- **Be careful not to overload the carrier especially with loose items as these may lodge in the wheels causing instability.**



Should any problems arise due to faulty manufacturing & materials under normal everyday use then it is covered by our comprehensive 12 months warrantee from the date of purchase (For terms and conditions please consult our Customer Services Department).

Specialised Orthotic Services Ltd cannot be held responsible for any unauthorised use or incorrect fitting whether resulting in personal injury or not and any damage caused to equipment not supplied with the tray.

The company has a policy of continual development and we reserve the right to alter specification and design without prior notification. All products manufactured are covered by design rights or patents. This document is the copyright of the company and the contents are its intellectual property.

If you require further copies of this handout or require further details relating to any of its content, then please contact us (Copyright applies):



Unit 127/128, Fauld Industrial Park, Tutbury, Staffordshire. DE13 9HS.

Tel: 01283 520400 Fax: 01283 520401

E-mail: enquiries@specialorthotic.com